

COMBOS



DIGI - COMBO 3 in 1



41 Pre Programmed

Size : 17x15.2x9.2 cm
 Net Wt. : 2.290 kg
 Shipping Wt. : 5 kg

Power Supply: 210 - 230 VAC 50Hz

Fuse Rating : 1A ±10%

Max. Power Consumption :

20W for DIGI - COMBO 3 in 1 & COMBO MST + TENS

55W for COMBO US + MST + TENS

Safety Standard IEC60601 - 1

Safetyclass Class I type BF (IEC60601 - 1)

SPECIFICATION

Ultra - Reiz Current

Frequency : 143 Hz
 Pulse width on time : 2m Sec.
 Pulse width off time : 5 m Sec.
 Output : 0- 160 V
 Pre Programmed : 11 pre programmed condition wise
 Patient safety : Output is 100% isolation
 Treatment Time : Digital Timer pre settable 30 minutes - 1 minute
 Display : Attractive Back light Bluish LCD
 Therapy Mode : G, IG, F and SF

Muscle Stimulator

Pulse duration : 0.01, 0.03, 0.1, 0.3, 1,3,10,30,300ms
 Pulse Rest Time : 1- 6 Sec. linearly adjustable
 Faradic Frequency : 50 Hz/ 0.7 m Sec.
 Surge period : 1- 6 Sec. linearly adjustable
 Surge Rest Time : 1- 6 Sec. linearly adjustable
 Output : 0 - 160 V

Tens

30 pre programmed condition wise
 Ref. DIGI TENS - Dual



COMBO US + MST + TENS



4 Patients can be treated at a time

Size : 24.5x28x10 cm
 Net Wt. : 4.5 kg
 Shipping Wt. : 7.5 kg



COMBO MST + TENS



Size : 15.3x16x7.7 cm
 Net Wt. : 2 kg
 Shipping Wt. : 3 kg

SPECIFICATION

Muscle Stimulator

Therapy mode : G, IG, F and SF
 Pulse duration : 1G Variable
 Visual Indicator : G, IG, F and SF
 Faradic Frequency : 50Hz/0.7 m Sec
 Surge Duration : Variable

Tens

Frequency : 2Hz to 100 Hz
 Pulse Width : 400 μ Sec
 Patient Safety : Output is 100% Isolation
 Therapy mode : continuous/burst
 Output : 0-100V
 Frequency Indicator : LED Indicator

Ultra Sound

Ref : Digi Sound - Pro or Digisound - Dual