

# ULTRA SOUND THERAPY

## PHYSIO SOUND - Plus



Size : 22.5x23.5x9.5 cm  
 Net Wt. : 2 kg  
 Shipping Wt. : 5 kg

### SPECIFICATION

Output Frequency : 1 MHz  
 Output Mode : Continuous & pulsed  
 Pulsed Mode : 1:2, 1:4, 1:8, 1: 16  
 Transducer : Normal Probe & Pointer Probe  
 Output power : Max 3.1 Watt/Cm<sup>2</sup> (Linearly Adjustable)  
 Patient Safety : Output is 100% isolated  
 Time : 1-60 Mts  
 Provision for combination therapy facility  
 Pointer Probe  
 ABS Cabinet)

Universal Power Supply: 110 - 230 VAC 50Hz / 60Hz  
 Fuse Rating : 3A ±10%  
 Max. Power Consumption : 35W  
 Safety Standard IEC60601 - 1  
 Safetyclass Class 1 type BF (IEC60601 - 1)

**Compact Size, Easy to carry for home visit**

## PHYSIO SOUND - Compact



Size : 16x16.5x7.5 cm  
 Net Wt. : 1.800 kg  
 Shipping Wt. : 4.800 kg

## DIGI SOUND - DUAL



**70 Pre Programmed Condition wise  
 Sperate Transducer for 1Mhz & 3 Mhz**

Size : 22.5x23.5x9.5 cm  
 Net Wt. : 2.500 kg  
 Shipping Wt. : 5.500 kg

### SPECIFICATION

Output mode : Continuous and pulsed  
 Pulse mode : 1:2, 1:4, 1:8, 1:16  
 Transducer : Separate Transducer for 1 MHz & 3 MHz  
 Indicator for Active Transducer : Light will glow in the active transducer  
 Parameter selection : Manual & Programmed  
 Intensity : Linear adjustment 3-1 w/cm<sup>2</sup>  
 Patient safety : Output is 100% isolation  
 Pre program : 70 pre programmed Condition Wise  
 for 1 MHz 40, for 3 MHz 30  
 Output display : Attractive Back Light Blueish LCD  
 Treatment Time : Digital timer Pre Settable 15 to 1 minute  
 Provision for Combination Therapy facility.  
 ABS Cabinet

## DIGI SOUND - Pro



**70 Pre Programmed Condition wise  
 1 MHz with two Transducer (pointer & Normal)**

Size : 22.5x23.5x9.5 cm  
 Net Wt. : 2.800 kg  
 Shipping Wt. : 5.800 kg