

# TENS



## TENS - 4 Channel



Size : 29.5x21.5x9 cm  
 Net Wt. : 3 kg  
 Shipping Wt. : 6 kg



## DIGI TENS - pro



30 Pre Programmed Condition wise

Size : 22.5x23.5x9.5 cm  
 Net Wt. : 2.690 kg  
 Shipping Wt. : 6.690 kg

### SPECIFICATION

Timer : Presetable 15 Mts  
 Display : Numeric LCD Display (DIGI TENS)  
 Frequency : 2 Hz to 250 Hz  
 Pulse width : 400 Micro Sec.  
 Therapy mode : All Channels continuous  
                   All Channels burst  
                   2 Channels continuous  
                   2 Channels burst  
 Output : Four Channel  
 Auto mode : Linear, Trapezoidal and Non Linear (DIGI TENS)  
                   (All the frequencies change one by one within 3 Sec.)  
 Parameter Selection: Manual & Programmed  
 Intensity : Linear Adjustment  
 Patient Safety : Output is 100% isolation  
 INTENSITY CAN BE SELECT SEPARATELY FOR CH1, CH2, CH3 and CH4  
 Power Supply : 210-230AC / 50Hz  
 Fuse Rating : 1A ±10%  
 Max. Power Consumption : 10W  
 Safety Standard IEC60601 - 1  
 Safetyclass Class 1 type BF (IEC60601 - 1)

### SPECIFICATION

Frequency : 2Hz to 150 Hz  
 Pulse width : 400 micro sec  
 Therapy Mode : Continuous and burst  
 Output : Two Channel  
 Timer (DIGI TENS) : 15 mts. in built timer.  
 Auto Mode : Linear, Trapezoidal & Non Linear (DIGI TENS)  
                   (all the frequencies change one by one with in 3 Sec.)  
 Intensity : linear adjustment  
 Patient Safety : Output is 100% isolation  
 Frequency Indicator : LED Indicator  
 INTENSITY CAN BE SELECT SEPARATELY FOR CH1 AND CH2  
 Power Supply : 9V battery or 9 V eliminator  
 Safety Standard IEC60601 - 1  
 Safetyclass Class 1 type BF (IEC60601 - 1)

## POCKET TENS - Dual



## DIGI TENS - Dual

