

# INTERFERENTIAL THERAPY



## PHYSIOMED - Plus



Size : 36x27x11 cm  
 Net Wt. : 3.650 kg  
 Shipping Wt. : 6.500 kg



## PHYSIOMED - Plus (Compact)



Size : 25.5x21.5x10 cm  
 Net Wt. : 2.350 kg  
 Shipping Wt. : 5.500 kg

### SPECIFICATION

Medium Frequency : 2 KHz, 4 KHz, 6 KHz, 8 KHz and 10 KHz  
 Therapy mode : 4p, 4pv, 2p  
 Vector : 4 types, 45°, 90° linear, 90° trapezoid, 90° non linear  
 Base : 0-150Hz  
 Spectrum : 1-100Hz  
 Output Display : Separate display for CH1 & CH2 Current, LED indicator for stimulation to the patient  
 Treatment time : Auto cutoff after 15 minutes  
 Spectrum changes : Bar Graph indicator  
 Spectrum programmer : Linear, Trapezoidal & Non linear  
 Isolation : Output 100 % isolated  
 Timer : Digital timer 15 minutes pre settable  
 Spectrum changes: Bar Graph indicator  
 TENS and Russian Current

### Compact Size Easy to Carry for Home Visit

Universal Power Supply: 110 - 275 VAC 50Hz / 60Hz  
 Fuse Rating : 1A ±10%  
 Max. Power Consumption : 20W  
 Safety Standard IEC60601 - 1  
 Safteyclass Class 1 type BF (IEC60601 - 1)



## DIGI IFT - Pro +



### SPECIFICATION

Medium Frequency : 2 KHz, 4 KHz, 6 KHz, 8 KHz and 10 KHz  
 Therapy mode : 4p, 4pv, 2p  
 Base : 1- 250 Hz  
 Sweep : 0- 250 Hz  
 Sweep Program : Linear, Trapezoidal, Saw tooth & Non linear  
 Patient Safety : Output is 100% Isolation  
 Therapy Mode : 4 P, 4 PV, 2P  
 Parameter Selection : Manual & Programmed  
 Display : Large Blueish Graphics LCD  
**RUSSION CURRENT**  
 Output : 2500 Hz Sine wave  
 ON time : 0 - 99 Sec.  
 OFF time : 0 - 99 Sec.  
 100 Pre Programmed Pain condition wise  
 100 User Memory with Patient name  
 Cable & Electrode Checking mode  
 ABS Cabinet

### IFT + Russian Current

100 Pre Programmed  
 100 User Memory  
 Compact Size  
 Easy to carry for home visit

Size : 22.5x23.5x9.5 cm  
 Net Wt. : 2.520 kg  
 Shipping Wt. : 4.500 kg