

VACUUM THERAPY

VACUUM - Plus



SPECIFICATION

- Vacuum Output : 0 to 600mm.hg/m Bar
- Vacuum Control : Linearly Adjustable.
- Mode : Continuous/Pulsed
- Pulsed Mode : Seven Different Pulses.
- Water Tap Full : Visual & Audio Indication.
- Input : Two RED, Two BLACK Socket (For Electrotherapy Combination Purpose)
- Power Supply : 200 - 230 VAC 50Hz
- Fuse Rating : 3A ± 10%
- Max. Power Consumption : 30W
- Safety Standard IEC60601 - 1
- Safetyclass Class II type B (IEC60601 - 1)

Size : 26x20.5x15 cm
 Net Wt. : 4 kg
 Shipping Wt. : 7 kg

MUSCLE STIMULATOR

PHYSIOSTIM-T

SPECIFICATION

- Therapy Mode : G, IG, F & SF
- Pulse duration : 0.01, 0.03, 0.1, 0.3, 1, 3, 10, 100, 300 and 1000m Sec.
- Pulse Rest time : 1-6 Sec.
- Visual Indicator : G, IG, F, & SF
- Faradic Frequency : 50Hz/0.7 m Sec.
- Surge Duration : 1-6 Sec.
- Surge Rest Time : 1-6 Sec.
- Power Supply : 200 - 230 VAC 50Hz
- Fuse Rating : 1A ± 10%
- Max. Power Consumption : 15W
- Safety Standard IEC60601 - 1
- Safetyclass Class 1 type BF (IEC60601 - 1)



PHYSIOSTIM-D



Size : 14x25x10 cm
 Net Wt. : 2.450 kg
 Shipping Wt. : 5 kg

Display : 3 digit Digital display meter available for diagnostic muscle stimulator