

PATIENT LIFT



Size : 50x36x69 INCH
 Net Wt. : 60 Kg
 Shipping Wt. : 80 Kg

SPECIFICATION

Base is 2" OD hollow square tube with powder coating surface.

The base is removable for easy storage and transportation

Can lift the patient weight upto 180 Kg Capacity with remote control

Four castors will provides easy mobility rear wheel brake will help to hold the system at any surface without movement.

Emergency to stop actuator movement

Operating voltage : built in rechargeable battery (24vDC)

Power Supply : 12 VDC

Fuse Rating : 1A \pm 10%

Max. Power Consumption : 5W

Safety Standard IEC60601 - 1

Safetyclass Class II type BF (IEC60601 - 1)

Electrical and Functional safety test certified

NEURO MUSCULAR CO ORDINATION THERAPY

SPECIFICATION

Separate displays for level and efficiency

L1 to L4 - Training Levels.

L5 to L7 - Testing levels.

Display shows efficiency in %

16 Touch pads with LEDs, (4x4)

Operation

Training Level.

On time L1-2Sec, L2-1.5 Sec, L3-1Sec, L4-0.5 Sec.

Each level with 3 Mts Running time.

Finally Display will shows efficiency in %

Testing Level.

On time & running time based on your speed.

Counting will stops when you fail to touch block.

Finally Display will shows efficiency in %

Power Supply : 110 - 230 VAC 50Hz / 60Hz Universal

Fuse Rating : 3A \pm 10%

Max. Power Consumption : 20W

Safety Standard IEC60601 - 1

Safetyclass Class 1 type BF (IEC60601 - 1)



Size Available :

1. 1x1.5 FEET : Childrens and Paralyzed Peoples

2. 4x4 FEET : Whole Body Exercises for Childrens

3. 8x8 FEET : Whole Body Exercises for Adult

Net Wt. : 4 kg Shipping Wt. : 6 kg

Net Wt. : 40 kg Shipping Wt. : 110 kg

Net Wt. : 60 kg Shipping Wt. : 120 kg