

# ELECTRONIC TRACTION



## PHYSIOTRAC - Plus



Size : 25x28.5x25.5 cm  
 Net Wt. : 9.300 kg  
 Shipping Wt. : 13 kg

### Imported AC Motor

Provide smooth operation without motor noise & jerk free traction.

### SPECIFICATION

Digital treatment timer  
 Separate Traction force for CERVICAL & LUMBAR.  
 Traction force for CERVICAL 4 Kg's to 15 Kg's (1Kg steps)  
 Traction force for LUMBAR 20 Kg's to 45Kg's  
 (2 Kg's steps) with Doublers up to 90 Kg's (optional)  
 LED Indicator for HOLD & REST time  
 Patient Safety Switch  
 Hold time 10,20,40,60 and 80 Sec.  
 Rest time 1, 5,10,15,20 Sec.  
 Static and Intermittent Mode

Power Supply: 210-230 VAC 50Hz  
 (Need 0.5 KVA Manual Stabilizer for Uniform Pulling)  
 Fuse Rating : 1A ±10%  
 Max. Power Consumption : 60W  
 Safety Standard IEC60601 - 1  
 Safetyclass Class II type B (IEC60601 - 1)

### SPECIFICATION

Display : Large LCD Display  
 Treatment Time : 1 to 99 mints  
 Rest Time : 1 to 99 sec  
 Hold Time : 1 to 99 sec  
 LED Indicators for Rest Time & Hold Time  
 LED Indicators for Static and Intermittent  
 50 User Programme  
 50 Pre Programme with Treatment Condition  
 Soft Touch Operation  
 Separate Traction force for CERVICAL & LUMBAR.  
 Traction force for CERVICAL 4 Kg's to 15 Kg's (1Kg steps)  
 Traction force for LUMBAR 20 Kg's to 45Kg's (2 Kg steps)  
 with Doubler upto 90 Kgs (optional)  
 Patient Safety Switch



## PHYSIOTRAC - Plus (Digital)



Size : 25x28.5x25.5 cm  
 Net Wt. : 9.300 kg  
 Shipping Wt. : 13 kg

## TRACTION TABLE



### SPECIFICATION

Total Length - 8 Feet, Height - 2½ Feet, Width - 2 Feet,  
 Bed Length - 6½ Feet, (Movable 3½ Feet)  
 Machine Fixing Stand Height Can be adjust  
 Net Wt. : 30 kg  
 Shipping Wt.: 40 kg