



LASER THERAPY (SCANNING and PORTABLE)

BMS LASER - Pro +

- Gives Guaranteed Result

Universal Power Supply : 110 - 275 VAC 50Hz / 60Hz
Fuse Rating : 1A $\pm$ 10%

Max. Power Consumption : 30W

Safety Standard IEC6061 - 1

Safetyclass Class 1 type BF (IEC6061 - 1)

Electrical and functional safety test certified



Size : 75x737x160cm
Net Wt. : 20 kg
Shipping Wt. : 50 kg

SPECIFICATION

Functions as both a Scanning Laser Therapy System and a Portable Laser Therapy Unit, offering flexible treatment options.

Laser arm adjustment: 0°–60° upward and 0°–60° downward, with optional remote-controlled positioning and 360° arm rotation for precise beam placement.

Both the Scanning Laser Control Unit and Portable Laser Therapy Unit feature premium mechanisms and a built-in accessory storage compartment for user convenience.

Digital control panel enables independent selection of laser beam pattern, scanning mode, and treatment area size, with a maximum scanning area of 2 ft along the both the X-axis and Y-axis.

Supports multiple beam patterns, including Spot, Circular, Elliptical, Horizontal, Vertical, and Combined Scanning Modes.

Automatic scanning ensures uniform laser energy distribution across the selected treatment area, promoting effective therapy and faster healing.