

# SHORT WAVE DIATHERMY



**SWD 500w**



**SWD 500 W  
Continuous / pulsed**

Size : 47x36x93 cm  
Net Wt. : 39 kg  
Shipping Wt. : 60 kg

## SPECIFICATION

Frequency : 27.12 MHz,  
Wave Length : 11.062 Meters  
Mode : Continuous/Pulsed  
Intensity : Adjustable in 5 steps  
Electrodes : One pair condenser pads with felt & cover  
Tuning : Manual tuning control knob to adjust the radio frequency for maximum output  
Timer : Digital timer pre settable 30 minutes at the end of the treatment announce the beep sound and automatic cut off output

Pulsed Mode  
Indicator provided for mains and RF  
Built in heavy duty cooling fan.  
Separate high frequency cable for disc electrodes with Arm  
Separate high frequency cable for pad electrodes.  
Heavy duty equipment mounted on castor to move place to place easily.  
Patient safety: output is 100% isolated  
Sinus pad, Disc electrode with Arms. (Optional)  
833 A Valve Used in 300 Watts & 500 Watts.



**SWD 300w**



**SWD 300 W  
Continuous / pulsed**

Size : 42x33x31 cm  
Net Wt. : 23 kg  
Shipping Wt. : 35 kg

Power Supply: 210 - 230 VAC 50Hz  
Fuse Rating : 5A ±10%  
Max. Power Consumption : 150W  
Safety Standard IEC60601 - 1  
Safetyclass Class I type BF (IEC60601 - 1)



# WAX BATH

## SPECIFICATION

Filament : "M" Type  
Thermostat : Deluxe Model  
Power Indicator : Available  
Thermostat Indicator : Available  
Safety fuse : Provided  
Thermostat control : 1 - 110° C  
Power Supply : 210 - 230 VAC 50Hz  
Fuse Rating : 5A ±10%  
Max. Power Consumption : 1200 W



**(Available in 1 kg 3 kgs 5 kgs 15 kgs & 30 kgs)**