



PARTNER FOR PHYSIO & REHABILITATION

# BMS

# Physiotherapy Equipments

- Designed for optimal performance

## UNDERWATER TREADMILL SYSTEM HYDROTHERAPY & FITNESS

Turn a one hour workout into 20 minutes, with better results and less strain on joints and muscles.

### QUICK FACTS

Overall Length	2500mm	98 inches
Overall Width	1150mm	45 inches
Overall Height	1800mm	71 inches
Treadmill Speed	0.5-12 Kph	

CE Approved Electricals


### EXERCISE CHAMBER

Inside Length	2000mm	79 inches
Inside Width	820mm	32 inches
Inside Height	1500mm	59 inches
Water Depth	0-1200mm	47 inches
Fill/Drain Time	2-3 Minutes	
Running Surface	1613mm	64 inches
Door Opening	673mm	26 inches
Weight	Empty: 411 Kg	Full: 1188 Kg.

### ELECTRICALS

208-240 VAC / 1 Phase/ 60 Hz/ 50 Amp Service

### AQUA REHAB TREADMILL



Aqua Rehab Treadmill

for multiple clinics and rehab centers. The exercise chamber gives additional privacy, and the larger track and adjustable console allow for customizable fitness and rehab needs. The Aqua Treadmill is an exceptional tool that utilizes the benefits of exercising in water offering low impact (buoyancy), with higher outcomes (resistance) and improved circulation (hydrostatic pressure)

### OPTIONAL FEATURES

- 10" Touchscreen Control Panel
- Programmable Treatments
- Foldable Deck for easy cleaning
- Resistance Jets
- Water Filtration and Temperature Control
- Can be kept outdoor or in hydro pool

**MEDICALLY APPROVED**

PAN INDIA SERVICING NETWORK

ON SITE TRAINING

24/7 SERVICE SUPPORT

## BHARAT MEDICAL SYSTEMS

AN ISO 9001 - 2008, ISO 13485-2012, CE, IEC 60601-1 Standard Certified

(Manufacturer, Exporter of Electro Therapy and Exercise Therapy Equipments)

M-1074, 7TH MAIN ROAD, THIRUVALLUVAR NAGAR, THIRUVANMIYUR, CHENNAI - 600 041, INDIA.

Phone No : +91 - 44 - 2451 6904 / 2451 6539

Whats app / Mobile No. : +91 - 94440 52532 / 98402 12531

90031 04470 / 99404 60676

## All Electro Therapy and Exercise Therapy Equipments Under One Roof



PARTNER FOR PHYSIO  
& REHABILITATION

# BMS

# Physiotherapy Equipments

- Designed for optimal performance

## THE COMPLETE BALANCE AND POSTURAL STABILITY ASSESSMENT SYSTEM



### Steps for Posture Balance Analysis

#### A. Baseline (Blind)

- Stand with both feet together for 60 seconds with eyes open.
- Stand with both feet together for 60 seconds with eyes closed.

#### B. Intervention (Balance Training)

Stand on the system for 3 minutes with eyes open while maintaining the posture.

- Maintain equal weight distribution on left and right feet.
- Maintain the balance of both the feet in anterior- posterior direction. Repeat the trial three times.
- Evaluation: Repeat both the steps of A again for the evaluation post intervention.



## **BHARAT MEDICAL SYSTEMS**

AN ISO 9001 - 2008, ISO 13485-2012, CE, IEC 60601-1 Standard Certified

(Manufacturer, Exporter of Electro Therapy and Exercise Therapy Equipments)

M-1074, 7TH MAIN ROAD, THIRUVALLUVAR NAGAR, THIRUVANMIYUR, CHENNAI - 600 041, INDIA.

Phone No : +91 - 44 - 2451 6904 / 2451 6539

Whats app / Mobile No. : +91 - 94440 52532 / 98402 12531

90031 04470 / 99404 60676

## All Electro Therapy and Exercise Therapy Equipments Under One Roof